Assembling instruction

* Part list:

Tart list.							
Part #	Module name	Part name	Inventor file name or product	thumbnail			
1		sensor house	sensor-house-reflection_v5.ipt				
2	sensor module	CMOS sensor	Basler ace acA3800-14uc Color USB 3.0 Camera				
3		sensor holder	sensor-holder-v3.ipt				
4		sensor slider	slave-slider-house.ipt				
5		linear ball bearing	LM8UU				
6	motor module	motor holder	master-slider-motor_v2.ipt				
7		motor	NEMA17				
8	mirror module of the "motor module	timing belt pulley wheel	GT2 Pully 20 teeth 5 mm	(
9		mirror holder of the "motor holder"	master-slider-mirror-v2.ipt				
10		belts	Mc master-Carr, High- Strength Ultra-Quiet Timing Belts				
11		frame holder	frame-holder-v3.ipt	a			
12		mirror framer holder	frame-holder-v3_MIR.ipt	6			
13		flanged ball bearing	uxcell F685ZZ				
14	sample holder module	incubator chamber holder	incubator-chamber-holder.ipt				
15		plate holder	plate-holder.ipt	\Diamond			


^{*} These instructions are for recreating the portable holographic imaging device used in "Stain-free, rapid, and quantitative viral plaque assay using deep learning and holography." Please feel free to contact liyuzhu@ucla.edu if you have any questions.

^{*} The dimensions of the 3D printed parts (.ipt) we provided should be further modified to accommodate the printing resolution of different 3D printers you will use.


16	light source module	light source holder part 1	LD_holder_p1_v2.ipt	
17		light source holder part 2	LD_holder_p2_vExdention.ipt	Ĵ
18		light source holder part 3	LD_holder_p3.ipt	
19		light source holder part 4	LD_holder_p4.ipt	J.
20		light source holder part 5	LD_holder_p5.ipt	

* Assembling steps:


Step 1: Finish the sensor module as follows.




Step 2: Finish the motor holder module as follows.




Step 3: Finish the mirror module of the "motor module."




Step 4: Connect the sensor module, the motor module, and its mirror module using two 8 mm rods (length: 250 mm). And then assemble the belt here.




Step 5: Add two 8 mm rods (length: 300 mm) onto the modules after "Step 3" as follows:




Step 6: Continue adding two frame holders (#11) and two mirror frame holders (#12) as follows:




Step 7: Continue to add four flanged ball bearings (#13) and two 5 mm rods (length: 250 mm) as follows:




Step 8: Continue adding five timing belt pulley wheels (#8) as follows:



Step 9: Continue assembling the rest of the belts and the second motor as follows. The 2D scanning stage is finished now.



Step 10: Continue to add the sample holder module.



Step 11: Continue to add the light source module. Now Finished!

